

# SOCKET GRAFTING

## MEMBRANE FREE CONE GRAFTS

If placing an implant, we recommend **Powerbone Cone**. If implant placement or timing is uncertain, we recommend **Biphasic Cone**.

### INCISION

Perform a mid-crest incision. Reflect/release sufficient to visualise site. A gently tissue stretch will assist closure.



### SITE PREP

Check for bone defects. Remove ALL degranulation tissue using a **Degranulation Burr**. A **Decortication** drill will assist healing.



### SOCKET CONE PREPARATION

Remove the socket cone from the packaging. Socket cone can be trimmed beforehand. For two rooted sites the socket cone can be partially split.



### SOCKET CONE APPLICATION

Apply the socket cone graft with fingers or forceps. Socket cone will absorb blood during application and swell to fill the socket. No membrane is needed.



### OPEN HEALING

A cross suture can be used to create tissue apposition. Apply **Ora Aid** dressing over the site. Alternatively, a sling type suture can be placed over the **Ora Aid** to prolongs coverage.



### IMPLANT PLACEMENT?

Place at 12 to 14 weeks following **Powerbone Socket Cone** (later if **Biphasic Socket Cone**). Check stability with **Penguin** & repeat ahead of loading.



Ivy Cottage, Bridgnorth Road,  
Stourton, DY7 5BH

M. +44 (0) 7800 990 131  
E. [ask@regen-store.co.uk](mailto:ask@regen-store.co.uk)  
W. [www.regen-store.co.uk](http://www.regen-store.co.uk)

Company registration number 10714405.



Incorporating Ariston Dental UK